



# BAHADIR STERILIZATION CONTAINERS- (STERRAD/Steris V-PRO Sterilization)

## INSTRUCTIONS OF USE

### Indications for Use:

Bahadir Sterilization Containers are reusable, metal, sterilization containers. They are designed for holding operating room instruments and/or textiles during low temperature hydrogen peroxide sterilization procedures and for maintaining sterility during storage and transport under proper hospital conditions. This container system is compatible for use with the following low temperature sterilizers and the cycles identified below:

\*STERIS V-PRO maX/maX2 - Lumen, Non-Lumen and Flexible Cycles

\*STERIS V-PRO s2/60 - Lumen, Non-Lumen and Flexible Cycles

\*STERRAD 100NX - Duo, Express, Standard and Flex Cycles

\*STERRAD NX - Standard and Advanced Cycles

### Bahadir USA Container Lumen Configurations

| Sterilization Cycle                     | Container Size     | Compatible Lumens   |
|---|--------------------|---|
| STERIS VPRO maX/maX2<br><br>Lumen Cycle | Full Size          | ≥0.77mm x ≤527mm (single, dual or triple channel)<br>≥0.8mm x ≤542mm (single, dual or triple channel)<br>≥0.48mm x ≤100mm (single, dual or triple channel) ≥1.3mm x ≤73mm (dead end lumen)<br>≥3mm x ≤298mm (rigid non-metallic lumen) ≥4mm x ≤424mm (rigid non-metallic lumen)       |
|   | Three-Quarter Size | ≥0.77mm x ≤527mm (single, dual or triple channel)<br>≥0.8mm x ≤542mm (single, dual or triple channel)<br>≥0.48mm x ≤100mm (single, dual or triple channel)<br>≥1.3mm x ≤73mm (dead end lumen)<br>≥3mm x ≤298mm (rigid non-metallic lumen)<br>≥4mm x ≤424mm (rigid non-metallic lumen) |
|   | Half Size          | ≥0.77mm x ≤527mm (single, dual or triple channel)<br>≥0.8mm x ≤542mm (single, dual or triple channel)<br>≥0.48mm x ≤100mm (single, dual or triple channel)<br>≥1.3mm x ≤73mm (dead end lumen)<br>≥3mm x ≤298mm (rigid non-metallic lumen)<br>≥4mm x ≤424mm (rigid non-metallic lumen) |

|   |                           |  |
|---|---------------------------|--|
| <p>STERIS VPRO maX/maX2</p> <p>Flexible Cycle</p> | <p>Full Size</p>          | <p>(2) flexible endoscopes with light cord (if not integral to endoscope) and mat with no additional load. The flexible endoscopes may contain Single or Dual channel lumens <math>\geq 1\text{mm} \times \leq 1050\text{mm}</math></p> <p>Or</p> <p>(1) flexible endoscope with light cord (if not integral to endoscope), endoscope accessories, mat and additional instruments. The flexible endoscope may contain Single or dual channel lumens that are <math>\geq 1\text{mm} \times \leq 1050\text{mm}</math>. Additional single, dual or triple channel stainless steel lumen device that is <math>\geq 0.48\text{mm} \times \leq 100\text{mm}</math> OR non-lumen device</p> |
|   | <p>Three-Quarter Size</p> | <p>(2) flexible endoscopes with light cord (if not integral to endoscope) and mat with no additional load. The flexible endoscopes may contain Single or Dual channel lumens <math>\geq 1\text{mm} \times \leq 1050\text{mm}</math></p> <p>Or</p> <p>(1) flexible endoscope with light cord (if not integral to endoscope), endoscope accessories, mat and additional instruments. The flexible endoscope may contain Single or dual channel lumens that are <math>\geq 1\text{mm} \times \leq 1050\text{mm}</math>. Additional single, dual or triple channel stainless steel lumen device that is <math>\geq 0.48\text{mm} \times \leq 100\text{mm}</math> OR non-lumen device</p> |
|   | <p>Half Size</p>          | <p>(2) flexible endoscopes with light cord (if not integral to endoscope) and mat with no additional load. The flexible endoscopes may contain Single or Dual channel lumens <math>\geq 1\text{mm} \times \leq 1050\text{mm}</math></p> <p>Or</p> <p>(1) flexible endoscope with light cord (if not integral to endoscope), endoscope accessories, mat and additional instruments. The flexible endoscope may contain Single or dual channel lumens that are <math>\geq 1\text{mm} \times \leq 1050\text{mm}</math>. Additional single, dual or triple channel stainless steel lumen device that is <math>\geq 0.48\text{mm} \times \leq 100\text{mm}</math> OR non-lumen device</p> |

|   |                    |  |
|---|--------------------|--|
| STERIS VPRO maX/maX2<br>Non-Lumen Cycle | Full Size          | -Non lumened instruments   |
|   | Three-Quarter Size | -Non lumened instruments   |
|   | Half Size          | -Non lumened instruments   |
| STERIS VPRO s2/60<br>Lumen cycle        | Full Size          | <p>≥0.77mm x ≤410mm (single or dual channel)</p> <p>≥1.8mm x ≤542mm (single or dual channel)</p> <p>≥1.2mm x ≤275mm (triple channel)</p> <p>≥1.8mm x ≤310mm (triple channel)</p> <p>≥2.8mm x ≤317mm (triple channel)</p>   |
|   | Three-Quarter Size | <p>≥0.77mm x ≤410mm (single or dual channel)</p> <p>≥1.8mm x ≤542mm (single or dual channel)</p> <p>≥1.2mm x ≤275mm (triple channel)</p> <p>≥1.8mm x ≤310mm (triple channel)</p> <p>≥2.8mm x ≤317mm (triple channel)</p>   |
|   | Half Size          | <p>≥0.77mm x ≤410mm (single or dual channel)</p> <p>≥1.8mm x ≤542mm (single or dual channel)</p> <p>≥1.2mm x ≤275mm (triple channel)</p> <p>≥1.8mm x ≤310mm (triple channel)</p> <p>≥2.8mm x ≤317mm (triple channel)</p>   |
| STERIS VPRO s2/60<br>Flexible Cycle     | Full Size          | <p>-Non-lumened instruments and instruments with diffusion restricted spaces (such as the hinged portion of forceps and scissors)</p> <p>-(1) surgical flexible endoscope or bronchoscope with light cord (if not integral to the endoscope). The flexible endoscope may be single or dual lumen device with lumens that are ≥1mm x ≤990mm</p> |
|   | Three-Quarter Size | <p>-Non-lumened instruments and instruments with diffusion restricted spaces (such as the hinged portion of forceps and scissors)</p> <p>-(1) surgical flexible endoscope or bronchoscope with light cord (if not integral to the endoscope). The flexible endoscope may be single or dual lumen device with lumens that are ≥1mm x ≤990mm</p> |
|   | Half Size          | <p>-Non-lumened instruments and instruments with diffusion restricted spaces (such as the hinged portion of forceps and scissors)</p> <p>-(1) surgical flexible endoscope or bronchoscope with light cord (if not integral to the endoscope). The</p>  |

|                                      |                    |  |
|--------------------------------------|--------------------|--|
|                                      |                    | flexible endoscope may be single or dual lumen device with lumens that are $\geq 1\text{mm}$ x $\leq 990\text{mm}$ |
|                                      |                    |  |
| STERIS VPRO s2/60<br>Non-Lumen Cycle | Full Size          | -Non lumened instruments   |
|                                      | Three-Quarter Size | -Non lumened instruments   |
|                                      | Half Size          | -Non lumened instruments   |
|                                      |                    |  |
| STERRAD NX<br>Standard Cycle         | Full Size          | (5) SS lumens 1mm I.D. x 150mm L<br>(5) SS lumens 2mm I.D. x 400mm L   |
|                                      | Three-Quarter Size | (10) SS lumens 1mm I.D. x 150mm L  |
|                                      | Half Size          | (10) SS lumens 1mm I.D. x 150mm L  |
|                                      |                    |  |
| STERRAD NX<br>Advanced Cycle         | Full Size          | (10) SS lumens 2mm I.D. x 500mm L  |
|                                      | Three-Quarter Size | (6) SS lumens 1mm I.D. x 350mm L   |
|                                      | Half Size          | (10) SS lumens 1mm I.D. x 200mm L  |
|                                      |                    |  |
| STERRAD 100NX<br>Standard Cycle      | Full Size          | (5) SS lumens 1mm I.D. x 500mm L   |
|                                      | Three-Quarter Size | (5) SS lumens 1mm I.D. x 350mm L   |
|                                      | Half Size          | (5) SS lumens 0.7mm I.D. x 200mm L   |
|                                      |                    |  |
| STERRAD 100NX<br>Express Cycle       | Full Size          | Non lumened SS instruments   |
|                                      | Three-Quarter Size | Non lumened SS instruments   |
|                                      | Half Size          | Non lumened SS instruments   |
|                                      |                    |  |
| STERRAD 100NX<br>Duo Cycle           | Full Size          | (1)PTFE tubing 1.5 mm I.D. x 850mm L   |
|                                      | Three-Quarter Size | (1)PTFE tubing 1.5 mm I.D. x 850mm L   |
|                                      | Half Size          | (1)PTFE tubing 1.5 mm I.D. x 850mm L   |
|                                      |                    |  |
| STERRAD 100 NX<br>Flexible Cycle     | Full Size          | (1) Scope lumen 1mm I.D. x 850mm L   |
|                                      | Three-Quarter Size | (1) Scope lumen 1mm I.D. x 850mm L   |
|                                      | Half Size          | (1)PTFE tubing 1mm I.D. x 850mm L  |

Bahadir USA Container Configurations

| Sterilization Cycle                        | Container Size            | Container Name | Total Loaded Container (lbs) |
|--|---------------------------|----------------|------------------------------|
| STERIS VPRO<br>maX/maX2<br>Lumen Cycle     | Full Size 4 1/8"          | Y111.10        | 19.65                        |
|  | Full Size 5 1/4"          | Y111.13        | 19.65                        |
|  | Full Size 5 7/8"          | Y111.15        | 19.65                        |
|  | Full Size 7 3/4"          | Y111.20        | 19.65                        |
|  | Full Size 10 1/4 "        | Y111.26        | 19.65                        |
|  | Full Size 7 3/4"          | Y111.70        | 19.65                        |
|  | Three-Quarter Size 4 1/8" | Y211.10        | 19.65                        |
|  | Three-Quarter Size 5 1/4" | Y211.13        | 19.65                        |
|  | Three-Quarter Size 5 7/8" | Y211.15        | 19.65                        |
|  | Three-Quarter Size 7 3/4" | Y211.20        | 19.65                        |
|  | Half Size 4 1/8"          | Y311.10        | 15                           |
|  | Half Size 5 1/4"          | Y311.13        | 15                           |
|  | Half Size 5 7/8"          | Y311.15        | 15                           |
|  | Half Size 7 3/4"          | Y311.20        | 15                           |
|  | Half Size 10 1/4"         | Y311.26        | 15                           |
| STERIS VPRO<br>maX/maX2<br>Flexible Cycle  | Full Size 4 1/8"          | Y111.10        | 21                           |
|  | Full Size 5 1/4"          | Y111.13        | 21                           |
|  | Full Size 5 7/8"          | Y111.15        | 21                           |
|  | Full Size 7 3/4"          | Y111.20        | 21                           |
|  | Full Size 10 1/4 "        | Y111.26        | 21                           |
|  | Full Size 7 3/4"          | Y111.70        | 21                           |
|  | Three-Quarter Size 4 1/8" | Y211.10        | 21                           |
|  | Three-Quarter Size 5 1/4" | Y211.13        | 21                           |
|  | Three-Quarter Size 5 7/8" | Y211.15        | 21                           |
|  | Three-Quarter Size 7 3/4" | Y211.20        | 21                           |
|  | Half Size 4 1/8"          | Y311.10        | 21                           |
|  | Half Size 5 1/4"          | Y311.13        | 21                           |
|  | Half Size 5 7/8"          | Y311.15        | 21                           |
|  | Half Size 7 3/4"          | Y311.20        | 21                           |
|  | Half Size 10 1/4"         | Y311.26        | 21                           |
| STERIS VPRO<br>maX/maX2<br>Non-Lumen Cycle | Full Size 4 1/8"          | Y111.10        | 26                           |
|  | Full Size 5 1/4"          | Y111.13        | 26                           |
|  | Full Size 5 7/8"          | Y111.15        | 26                           |
|  | Full Size 7 3/4"          | Y111.20        | 26                           |
|  | Full Size 10 1/4 "        | Y111.26        | 26                           |
|  | Full Size 7 3/4"          | Y111.70        | 26                           |
|  | Three-Quarter Size 4 1/8" | Y211.10        | 26                           |
|  | Three-Quarter Size 5 1/4" | Y211.13        | 26                           |
|  | Three-Quarter Size 5 7/8" | Y211.15        | 26                           |
|  | Three-Quarter Size 7 3/4" | Y211.20        | 26                           |
|  | Half Size 4 1/8"          | Y311.10        | 19                           |
|  | Half Size 5 1/4"          | Y311.13        | 19                           |
|  | Half Size 5 7/8"          | Y311.15        | 19                           |
|  | Half Size 7 3/4"          | Y311.20        | 19                           |
|  | Half Size 10 1/4"         | Y311.26        | 19                           |
|  | Full Size 4 1/8"          | Y111.10        | 11                           |
|  | Full Size 5 1/4"          | Y111.13        | 11                           |

|   |                           |         |      |
|---|---------------------------|---------|------|
| STERIS VPRO<br>S2/60<br>Lumen Cycle     | Full Size 5 7/8"          | Y111.15 | 11   |
|   | Full Size 7 3/4"          | Y111.20 | 11   |
|   | Full Size 10 1/4 "        | Y111.26 | 11   |
|   | Full Size 7 3/4"          | Y111.70 | 11   |
|   | Three-Quarter Size 4 1/8" | Y211.10 | 11   |
|   | Three-Quarter Size 5 1/4" | Y211.13 | 11   |
|   | Three-Quarter Size 5 7/8" | Y211.15 | 11   |
|   | Three-Quarter Size 7 3/4" | Y211.20 | 11   |
|   | Half Size 4 1/8"          | Y311.10 | 11   |
|   | Half Size 5 1/4"          | Y311.13 | 11   |
|   | Half Size 5 7/8"          | Y311.15 | 11   |
|   | Half Size 7 3/4"          | Y311.20 | 11   |
|   | Half Size 10 1/4"         | Y311.26 | 11   |
|   |                           |         |      |
| STERIS VPRO<br>S2/60<br>Flexible Cycle  | Full Size 4 1/8"          | Y111.10 | 11   |
|   | Full Size 5 1/4"          | Y111.13 | 11   |
|   | Full Size 5 7/8"          | Y111.15 | 11   |
|   | Full Size 7 3/4"          | Y111.20 | 11   |
|   | Full Size 10 1/4 "        | Y111.26 | 11   |
|   | Full Size 7 3/4"          | Y111.70 | 11   |
|   | Three-Quarter Size 4 1/8" | Y211.10 | 11   |
|   | Three-Quarter Size 5 1/4" | Y211.13 | 11   |
|   | Three-Quarter Size 5 7/8" | Y211.15 | 11   |
|   | Three-Quarter Size 7 3/4" | Y211.20 | 11   |
|   | Half Size 4 1/8"          | Y311.10 | 11   |
|   | Half Size 5 1/4"          | Y311.13 | 11   |
|   | Half Size 5 7/8"          | Y311.15 | 11   |
|   | Half Size 7 3/4"          | Y311.20 | 11   |
| Half Size 10 1/4"                       | Y311.26                   | 11      |      |
|   |                           |         |      |
| STERIS VPRO<br>S2/60<br>Non-Lumen Cycle | Full Size 4 1/8"          | Y111.10 | 19   |
|   | Full Size 5 1/4"          | Y111.13 | 19   |
|   | Full Size 5 7/8"          | Y111.15 | 19   |
|   | Full Size 7 3/4"          | Y111.20 | 19   |
|   | Full Size 10 1/4 "        | Y111.26 | 19   |
|   | Full Size 7 3/4"          | Y111.70 | 19   |
|   | Three-Quarter Size 4 1/8" | Y211.10 | 19   |
|   | Three-Quarter Size 5 1/4" | Y211.13 | 19   |
|   | Three-Quarter Size 5 7/8" | Y211.15 | 19   |
|   | Three-Quarter Size 7 3/4" | Y211.20 | 19   |
|   | Half Size 4 1/8"          | Y311.10 | 12.5 |
|   | Half Size 5 1/4"          | Y311.13 | 12.5 |
|   | Half Size 5 7/8"          | Y311.15 | 12.5 |
|   | Half Size 7 3/4"          | Y311.20 | 12.5 |
| Half Size 10 1/4"                       | Y311.26                   | 12.5    |      |
|   |                           |         |      |
| STERRAD NX                              | Full Size 4 1/8"          | Y111.10 | 10.7 |
|   | Full Size 5 1/4"          | Y111.13 | 10.7 |
|   | Full Size 5 7/8"          | Y111.15 | 10.7 |
|   | Full Size 7 3/4"          | Y111.20 | 10.7 |
|   | Full Size 10 1/4 "        | Y111.26 | 10.7 |
|   | Full Size 7 3/4"          | Y111.70 | 10.7 |
|   | Three-Quarter Size 4 1/8" | Y211.10 | 10.7 |
|   | Three-Quarter Size 5 1/4" | Y211.13 | 10.7 |
|   | Three-Quarter Size 5 7/8" | Y211.15 | 10.7 |

|                                |                                 |           |         |         |      |
|--------------------------------|---------------------------------|-----------|---------|---------|------|
| Standard Cycle                 | Three-Quarter Size              | 5 7/8"    | Y211.15 | 10.7    |      |
|                                | Three-Quarter Size              | 7 3/4"    | Y211.20 | 10.7    |      |
|                                | Half Size                       | 4 1/8"    | Y311.10 | 10.7    |      |
|                                | Half Size                       | 5 1/4"    | Y311.13 | 10.7    |      |
|                                | Half Size                       | 5 7/8"    | Y311.15 | 10.7    |      |
|                                | Half Size                       | 7 3/4"    | Y311.20 | 10.7    |      |
|                                | Half Size                       | 10 1/4"   | Y311.26 | 10.7    |      |
|                                |                                 |           |         |         |      |
| STERRAD NX<br>Advanced Cycle   | Full Size                       | 4 1/8"    | Y111.10 | 10.7    |      |
|                                | Full Size                       | 5 1/4"    | Y111.13 | 10.7    |      |
|                                | Full Size                       | 5 7/8"    | Y111.15 | 10.7    |      |
|                                | Full Size                       | 7 3/4"    | Y111.20 | 10.7    |      |
|                                | Full Size                       | 10 1/4"   | Y111.26 | 10.7    |      |
|                                | Full Size                       | 7 3/4"    | Y111.70 | 10.7    |      |
|                                | Three-Quarter Size              | 4 1/8"    | Y211.10 | 13.85   |      |
|                                | Three-Quarter Size              | 5 1/4"    | Y211.13 | 13.85   |      |
|                                | Three-Quarter Size              | 5 7/8"    | Y211.15 | 13.85   |      |
|                                | Three-Quarter Size              | 7 3/4"    | Y211.20 | 13.85   |      |
|                                | Half Size                       | 4 1/8"    | Y311.10 | 10.7    |      |
|                                | Half Size                       | 5 1/4"    | Y311.13 | 10.7    |      |
|                                | Half Size                       | 5 7/8"    | Y311.15 | 10.7    |      |
|                                | Half Size                       | 7 3/4"    | Y311.20 | 10.7    |      |
|                                | Half Size                       | 10 1/4"   | Y311.26 | 10.7    |      |
|                                |                                 |           |         |         |      |
|                                | STERRAD 100NX<br>Standard Cycle | Full Size | 4 1/8"  | Y111.10 | 21.4 |
| Full Size                      |                                 | 5 1/4"    | Y111.13 | 21.4    |      |
| Full Size                      |                                 | 5 7/8"    | Y111.15 | 21.4    |      |
| Full Size                      |                                 | 7 3/4"    | Y111.20 | 21.4    |      |
| Full Size                      |                                 | 10 1/4"   | Y111.26 | 21.4    |      |
| Full Size                      |                                 | 7 3/4"    | Y111.70 | 21.4    |      |
| Three-Quarter Size             |                                 | 4 1/8"    | Y211.10 | 13.85   |      |
| Three-Quarter Size             |                                 | 5 1/4"    | Y211.13 | 13.85   |      |
| Three-Quarter Size             |                                 | 5 7/8"    | Y211.15 | 13.85   |      |
| Three-Quarter Size             |                                 | 7 3/4"    | Y211.20 | 13.85   |      |
| Half Size                      |                                 | 4 1/8"    | Y311.10 | 13.85   |      |
| Half Size                      |                                 | 5 1/4"    | Y311.13 | 13.85   |      |
| Half Size                      |                                 | 5 7/8"    | Y311.15 | 13.85   |      |
| Half Size                      |                                 | 7 3/4"    | Y311.20 | 13.85   |      |
| Half Size                      |                                 | 10 1/4"   | Y311.26 | 13.85   |      |
|                                |                                 |           |         |         |      |
| STERRAD 100NX<br>Express Cycle |                                 | Full Size | 4 1/8"  | Y111.10 | 21.4 |
|                                | Full Size                       | 5 1/4"    | Y111.13 | 21.4    |      |
|                                | Full Size                       | 5 7/8"    | Y111.15 | 21.4    |      |
|                                | Full Size                       | 7 3/4"    | Y111.20 | 21.4    |      |
|                                | Full Size                       | 10 1/4"   | Y111.26 | 21.4    |      |
|                                | Full Size                       | 7 3/4"    | Y111.70 | 21.4    |      |
|                                | Three-Quarter Size              | 4 1/8"    | Y211.10 | 13.85   |      |
|                                | Three-Quarter Size              | 5 1/4"    | Y211.13 | 13.85   |      |
|                                | Three-Quarter Size              | 5 7/8"    | Y211.15 | 13.85   |      |
|                                | Three-Quarter Size              | 7 3/4"    | Y211.20 | 13.85   |      |
|                                | Half Size                       | 4 1/8"    | Y311.10 | 13.85   |      |
|                                | Half Size                       | 5 1/4"    | Y311.13 | 13.85   |      |
|                                | Half Size                       | 5 7/8"    | Y311.15 | 13.85   |      |
|                                | Half Size                       | 7 3/4"    | Y311.20 | 13.85   |      |

|                                     |                         |         |       |
|-------------------------------------|-------------------------|---------|-------|
|                                     | Half Size 10 ¼"         | Y311.26 | 13.85 |
|                                     |                         |         |       |
| STERRAD 100NX<br><br>DUO Cycle      | Full Size 4 ⅞"          | Y111.10 | 13.2  |
|                                     | Full Size 5 ¼"          | Y111.13 | 13.2  |
|                                     | Full Size 5 ⅞"          | Y111.15 | 13.2  |
|                                     | Full Size 7 ¾"          | Y111.20 | 13.2  |
|                                     | Full Size 10 ¼ "        | Y111.26 | 13.2  |
|                                     | Full Size 7 ¾"          | Y111.70 | 13.2  |
|                                     | Three-Quarter Size 4 ⅞" | Y211.10 | 13.2  |
|                                     | Three-Quarter Size 5 ¼" | Y211.13 | 13.2  |
|                                     | Three-Quarter Size 5 ⅞" | Y211.15 | 13.2  |
|                                     | Three-Quarter Size 7 ¾" | Y211.20 | 13.2  |
|                                     | Half Size 4 ⅞"          | Y311.10 | 7.2   |
|                                     | Half Size 5 ¼"          | Y311.13 | 7.2   |
|                                     | Half Size 5 ⅞"          | Y311.15 | 7.2   |
|                                     | Half Size 7 ¾"          | Y311.20 | 7.2   |
|                                     | Half Size 10 ¼"         | Y311.26 | 7.2   |
|                                     |                         |         |       |
| STERRAD 100NX<br><br>Flexible Cycle | Full Size 4 ⅞"          | Y111.10 | 12.59 |
|                                     | Full Size 5 ¼"          | Y111.13 | 12.59 |
|                                     | Full Size 5 ⅞"          | Y111.15 | 12.59 |
|                                     | Full Size 7 ¾"          | Y111.20 | 12.59 |
|                                     | Full Size 10 ¼ "        | Y111.26 | 12.59 |
|                                     | Full Size 7 ¾"          | Y111.70 | 12.59 |
|                                     | Three-Quarter Size 4 ⅞" | Y211.10 | 13.85 |
|                                     | Three-Quarter Size 5 ¼" | Y211.13 | 13.85 |
|                                     | Three-Quarter Size 5 ⅞" | Y211.15 | 13.85 |
|                                     | Three-Quarter Size 7 ¾" | Y211.20 | 13.85 |
|                                     | Half Size 4 ⅞"          | Y311.10 | 13.85 |
|                                     | Half Size 5 ¼"          | Y311.13 | 13.85 |
|                                     | Half Size 5 ⅞"          | Y311.15 | 13.85 |
|                                     | Half Size 7 ¾"          | Y311.20 | 13.85 |
|                                     | Half Size 10 ¼"         | Y311.26 | 13.85 |

Bahadir USA containers may be used with Bahadir USA accessories according to the table below:

| Cycle  | Baskets including holding pins, holding clamps, silicone holders, partition sheets, tamper evident locks | Silicone Mats |
|--|--|---------------|
| Steris VPRO 60/s2<br>Lumen, Non-Lumen, Flex    | YES  | YES           |
| Steris VPRO maX/maX2<br>Lumen, Non-Lumen, Flex | YES  | YES           |
| STERRAD NX – Standard                          | YES  | NO            |
| STERRAD NX – Advanced                          | YES  | NO            |
| STERRAD 100NX – Standard                       | YES  | NO            |
| STERRAD 100NX – Duo                            | YES  | NO            |
| STERRAD 100NX – Flex                           | YES  | NO            |
| STERRAD 100NX - Express                        | YES  | YES           |



## **DIRECTIONS FOR USE**

This user manual describes important instructions on the proper usage and maintenance of Bahadir Sterilization Containers, and possible hazards that could result from failure to observe the instructions.

- ❖ In order to prevent damage to the containers during usage it is important that the containers be handled by fully trained and qualified personnel in all areas of sterilization containers, hospital hygiene and sterilization technology.
- ❖ Thoroughly clean prior to first use. Items are shipped Non-Sterile

### **PRE-CLEANING**

- ✓ Remove lid and basket with instruments from the container.
- ✓ Separate the safety lid from the container lid, (if using safety lid)
- ✓ Remove filter retainer plate(s). (If bottom perforated, also remove retainer plates)
- ✓ Remove indicators and disposable security seals and disposable paper filter.

### **CLEANING AND DISINFECTION**

Requirements according DIN 58953/9; Users have to specify by means of disinfection and cleaning plan, when and how the sterilization containers have to be cleaned and/or disinfected. Containers that are used for waste disposal have to be cleaned and disinfected each time after use.



Do not use metal brushes or any abrasive cleaning material. Never use bleach or other corrosive chemicals to clean Bahadir Sterilization Containers. Abrasive products may cause chemical and physical corrosions that will void all manufacturer warranties.

For more information concerning the inspection and cleaning of closed sterilization containers, please see Sections 5.2.3 and 7.5.9 of ANSI/AAMI ST79 Standard, if available, or contact Bahadir USA at [info@bahadirusa.com](mailto:info@bahadirusa.com) or call 856-517-3080

### **MANUAL CLEANING**

See Table 1 for Manual Washing Parameters

1. Remove all remaining external process indicators and disposable locks
2. Remove lid from bottom of container
3. Remove the basket and any instruments from the container
4. Single-use Filter
  - a. Remove the Retention plate(s)
  - b. Remove filter(s) and discard (if present)
5. Rinse visible debris from all container components for at least 2 minutes under cold tap water
  - a. If debris remains, repeat rinsing
6. Prepare a bath of cleaning solution using STERIS® Prolystica® 2x Concentrate Neutral Cleaner according to manufacturer's instructions.
7. Fully immerse the container components in cleaning solution and soak for minimum of 5 minutes.
8. Use a soft bristled brush or sponge and scrub each component for approximately 5 minutes
9. Rinse all the parts under cold deionized water for 2 minutes without leaving any stain or residue on them.
10. Thoroughly dry with soft lint-free cloth or compressed filtered air and store in clean and dry environment.

**Table 1: Manual Washing Parameters:**

| <b>Cleaning Process Steps</b>  | <b>Instructions</b>  |
|--------------------------------|--|
| Prepare container for cleaning | Remove any indicators and locks, remove basket and any instruments, separate base from lid, remove filter retention plates, discard filters if present |
| Rinse                          | Using cold tap water, rinse all visible debris from components for at least 2 minutes.   |
| Cleaning Solution              | STERIS® Prolystica® 2x Neutral Cleaner, prepare bath according to manufacturer directions  |
| Soak                           | Fully immerse components for minimum for 5 minutes   |
| Cleaning/Brushing              | Scrub with soft sponge or brush to clean components for approximately 5 minutes  |
| Final Rinse                    | Thoroughly rinse for 2 minutes under cold deionized water to remove all detergent residues   |
| Dry                            | Dry using a lint free cloth or filtered compressed air   |
| Inspect                        | After cleaning visually inspect and repeat the cleaning process if a visually clean endpoint has not been achieved                                     |

**MECHANICAL CLEANING**

*Mechanical cleaning of the containers is preferred to manual cleaning provided that the washing machine is designed for cleaning sterilization containers. This applies in particular to ensure secure placement in the washing baskets and arrangement of the spray jets and arms. Washer must also have a special washing program for aluminum containers. See Table 2 for Mechanical Washing Parameters*



Do not use any cleaning solutions that contain soda or caustic acid. Do not use additional acidic neutralizers.

1. Remove all remaining external process indicators and disposable locks
2. Remove lid from bottom of container
3. Remove the basket and any instruments from the container
4. Single-use Filter
  - a. Remove the Retention plate(s)
  - b. Remove filter(s) and discard (if present)
  - c. Rinse visible debris from the retention plates
5. Rinse visible debris from all container components
  - a. Demineralized water is recommended for the final rinse and steam generation to avoid discoloration or damage resulting from minerals found in utility water
  - b. Make sure no residues from the Pre-cleaning process remain on the lid or bottom
  - c. If debris remains, repeat rinsing
6. Place components on washer rack facing down to avoid water collection
  - a. Retention plates should be placed in a separate washing basket
  - b. Use a neutral pH enzymatic cleaning solution in accordance with the manufacturer’s instructions
7. After mechanical cleaning cycle thoroughly dry (either with a soft, dry cloth or air dry) all components, including retention plate(s) before proceeding to assembly

**Table 2: Mechanical Washer Cycle Parameters:**

| Treatment   | Time (mm:ss) | Temperature          | Cleaning Solution                                     |
|-------------|--------------|----------------------|---|
| Pre-Rinse   | 01:30        | Cold Tap             | N/A   |
| Wash        | 05:00        | 55°C                 | STERIS® Prolystica® 2X<br>Concentrate Neutral Cleaner |
| Final Rinse | 01:45        | Deionized Water 90°C | N/A   |
| Dry         | 07:30        | 90°C                 | N/A   |

**ASSEMBLY FOR USE**

Bahadir Sterilization Containers must be checked visually before each usage. Inspect the container bottoms, container lids and the surfaces where seals sit must be free of dents, visible deformations and anything else that might appear out of the ordinary.

Users are responsible for pre-inspection prior to assembly for use of containment devices to verify that the gaskets in the lid and filter retention plates are in good working order, the plates properly seal over the vented area of the container lid and base, and that no rust or corrosion are present on any surfaces. Furthermore, the gasket should be intact and the edge of the container base nest properly within the container lid. Dents and deformations can affect the seal of the container compromising the integrity of the containment device. If any of the above issues exist, contact manufacturer.



Do not use spray oil or solvents on the lid seals. The seal in the inner lid must be completely inserted and undamaged. If there is any kind of damage detected, lids should not be used.

**PRE-ASSEMBLY INSPECTION AND PREPARATION:**


- ✓ Select surgical instrument set for sterilization.
- ✓ Be sure all the components are completely dry.
- ✓ Select appropriate size baskets and containers.

Bahadir Sterilization Containers are designed to be used with single-use Polypropylene filters validated for use with the containers (BAHADIR USA item number Y801.01).

**FILTER ASSEMBLY:**

1. Place one sheet of filter over each perforated section on the inside of Bahadir Sterilization Containers. (Make sure it covers all the perforation holes properly). User must ensure that only the approved polypropylene filter is used.
2. Make sure that filters and filter retainers are placed properly. ‘Click’ sound should be heard while placing filter retainers by pressing on them which indicates locking is realized.

\*Ensure the filter retention plates are seated properly and securely locked tight in position



#### **CAUTION:**

Filters and filter retainers must show no visual deformations. These parts must be checked visually and for their functionality before usage.

#### **INSTRUMENT AND CONTAINER ASSEMBLY**

Place thoroughly cleaned and dried instruments into the instrument basket(s) according to established hospital procedures. Complex instruments, such as air-powered instruments have to be prepared and sterilized according to the instrument manufacturer's instructions. Small baskets, trays and other types of accessories, especially with covers and lids, should only be used with the sterilization container which has been specifically designed for that purpose. Place instrument basket(s) into the prepared container bottom. Place assembled lid onto the container bottom, aligning handles on bottom with latches on lid. Simultaneously close both locking latches on the lid.

**SPECIAL Note:** Devices placed within the container are not to be separately enclosed in sterile barriers because inadvertent use of such secondary packaging by users has the potential to cause failures in the recommended sterilization cycles and raises the risk of infections for patients due to reduced sterility assurance levels.

**CAUTION:** Leave one inch of free space between the instruments and the inside of the container lid for effective processing.

#### **PROCESSING ASSEMBLY**

1. Select the appropriate indicator card and insert into the holding brackets on the outside of the container.
2. Insert the security seals (Y819.00).
3. Secure and Lock.

Bahadir Sterilization Containers are protected by disposable plastic security seals which once attached can be opened only by breaking which prevents inadvertent opening of the containers, and ensures that it is evident whether or not the sterility of the contents has been compromised. Use of internal and external indicators should be in accordance with in-house protocol that is determined by the user.

#### **LOADING THE STERILIZER:**

Bahadir Sterilization Containers should be placed horizontally in the sterilizer for effective sterilization and drying.

#### **STACKING:**

Stacking NOT permitted in the chamber during the sterilization

#### **PROCESSING**

The loading instructions of the sterilizer manufacturer should be observed

**CAUTION:** In order to avoid post-sterilization against accidents Containers need to remain on the container cart, in draft free area until cool enough to handle

Sterilization Cycle Parameters and Configurations for the Bahadir Container System by Modality – See Indications for Use at beginning of manual.

#### **SHELF LIFE**

In protected storage (temperature, humidity, air filtration etc. controlled hospital storage room conditions) containers can store sterilized medical devices for up to 180 days, provided the integrity of the container is not compromised (Tamper-evident security seal (Y819.00) remains intact). Facilities should follow industry guidelines and establish processes and procedures related to shelf-life and storage.

**Note:** Bahadir Sterilization Containers can be stacked during storage, one on top of another. Be sure that heavier container is placed at the bottom.

#### **BAHADIR USA CONTAINER SYSTEM WARRANTY**

BAHADIR USA hereby warrants all products carrying the registered trademark and logo of BAHADIR USA purchased directly from BAHADIR USA or an authorized dealer in accordance with the following terms and conditions: The warranty period of Bahadir Sterilization Containers is two years from the date of purchase. Bahadir Sterilization Containers are guaranteed to be free from defects in material and workmanship when used for its intended purpose. Containers that prove defective in material or workmanship will either be repaired/replaced free of charge during the warranty period. Containers that are not purchased from BAHADIR USA or an authorized source may not meet BAHADIR USA quality standards.

**Repair or Replacement:** Under warranty, BAHADIR USA will repair or replace any parts and determine the nature of any defect, the necessity and manner of repair or replacement, all other matters pertaining to the conditions of the Bahadir Sterilization Container Systems. All devices returned for repair or replacement are assessed on an individualized basis, and may result in either repair or replacement upon thorough inspection and evaluation with the safety and integrity of the device as our foremost concern. Repair will be limited to components which would not affect the performance of the containers with regards to sterilization, such as handles. Containers which have damage to the frame or seal areas which could affect the integrity of the system will be replaced.

#### **Exclusions**

This warranty shall not apply to any condition(s) or damage resulting from misuse, improper cleaning, improper handling, negligence, improper opening or unauthorized repair work including but not limited to: Incidents of abuse such as denting of Bahadir Sterilization Container due to dropping or other instances of mechanically applied pressure. Damage determined to be related to caustic or abrasive cleaning agents, items modified by the customer, items modified or customized by BAHADIR USA at the request of the customer, damage from fire, flood and other occurrences not under the control of BAHADIR USA.

Please direct your questions regarding this warranty to:

**BAHADIR USA**

431 South Pennsville Auburn Road

Carney's Point, NJ 08069

Tel: (856) 514-3080 Fax: (856) 514-3081

[www.bahadirusa.com](http://www.bahadirusa.com)

#### **Label Symbols Key:**



Product is non-sterile



Consult Instructions for Use



Caution



**Manufacturer:**

**BAHADIR USA LLC**

**431 South Pennsville Auburn Road**

**Carney's Point, NJ 08069**

**Tel: (856) 514-3080 Fax: (856) 514-3081**

**[www.bahadirusa.com](http://www.bahadirusa.com)**

INS00 Rev A